

Dental Mirrors



**Dental Mirror
Di-0001**

50mm Diameter 25D Value Range Mirror
New A new design equine dental mirror with silicon coated framed head and T handle. The shaft of the mirror is made up of stainless steel. Due to the combination of various materials this mirror is valued at a very competitive price.



**Dental Mirror
Di-0002**

75mm Diameter 25D Inclined Dental Mirror
Dental mirror for examination purposes of equine mouth.



**Dental Mirror
Di-0003**

50mm Diameter 25 Degree Inclined Dental mirror for examination purposes of equine mouth.



**Dental Mirror
Di-0004**


50mm Diameter Mirror (Straight) Dental mirror for examination purposes of equine mouth.



**Dental Mirror
Di-0005**

75mm Diameter Straight Head Mirror Dental mirror for examination purposes of equine mouth.

Dental Picks



**Dental Pick
Di-0006**

Lingual Buccal Long Head Dental Pick Its lingual and buccal application makes it ideal



**Dental Pick
Di-0007**

Long Head Dental Pick Set of 3



**Dental Pick
Di-0008**

(Anterior Posterior) Short Head Dental Pick For anterior and posterior use when elevating tissues around molar teeth, usually



**Dental Pick
Di-0009**



**Dental Pick
Di-0010**



**Dental Pick
Di-0011**

Dentistry Kit Packages



Dentistry Kit Package

Di-0012

Bronze Kit For Dentistry Kit includes Hausman speculum with biothene straps Float SH1 - 55cm premolar and molar float with 13 degree angled head Float SH5 - 55cm straight shaft float for cheek teeth Float 26 - 55cm upturned 15 degree, back molar float 3 corresponding blades for kit floats 1 x water proof wallet for floats and blades.



Dentistry Kit Package

Di-0013

Silver package for Equine Dentistry Kit includes SE speculum with biothene straps, a dental mirror, a 100ml syringe, a diamond S file and prime head light. Float SH1 - 55cm premolar and molar float with 13 degree angled head Float SH5 - 55cm straight shaft float for cheek teeth Float 26 - 55cm upturned 15 degree, back molar float 3 corresponding blades for the kit floats 1 x water proof wallet for floats and blades.



Dentistry Kit Package

Di-0014

Gold Package for Dentistry Kit includes Millennium speculum with extra wide bite plates, 100ml drenching syringe, bucket with brush, a dental mirror, 9.5" wolf tooth forceps ,a diamond S file and prime head light. Float SH1 - 55cm premolar and molar float with 13 degree angled head Float SH2 - 44cm premolar float with 15 degree angled head Float SH3 - 44cm offset flat shaft float Float SH5 - 55cm straight shaft float for cheek teeth Float 26 - 55cm upturned 15 degree, back molar float 5 corresponding blades for the kit floats 1 x water proof wallet for floats and blades.

Equine Mouth Speculums



**Equine Mouth Speculum
Di-0015**

Rubber Fitted Extra Wide Bite Plates Special rubber layer insertion provides protection against the chipping and slippage of the incisors. The extra wide arms keep the speculum away from the equine's cheek making more working space for the technicians to use various angles to access the horse's mouth. For use with the Millennium speculum, the SE Speculum or the Hausmans.



**Equine Mouth Speculum
Di-0016**

Large Extra Wide Bite Plates Made from premium stainless steel, these bite plates are flat and very deep, making them great for mis-aligned incisors. The extra wide arms keep the speculum away from the equine's cheek making more working space for the technicians to use various angles to access the horse's mouth. For use with Millennium speculum, the SE speculum or the hausman's.



**Equine Mouth Speculum
Di-0017**

Extra Wide Bite Plates This set of plates due to their semi circle shape prevent the incisors to come off the bite plate. The extra wide arms keep the speculum away from the equine's cheek making more working space for the technicians to use various angles to access the horse's mouth. For use with the Millennium speculum, the SE Speculum or the Hausmans.



**Equine Mouth Speculum
Di-0018**

Made up of polished forged stainless steel. Suitable for both Hausman Speculum and Millennium Speculum.



**Equine Mouth Speculum
Di-0019**

Incisor Bite Plates Specially designed silicon coated plates are suitable for the horses with damaged incisors. For use with the



**Equine Mouth Speculum
Di-0020**

Bite Plates For Pony Speculum Smaller size classical bite plates are specially designed for ponies or miniature breeds. Appropriate for all types of speculums.

Equine Mouth Speculums



**Equine Mouth Speculum
Di-0021**

SE Speculum with Biothene Straps 6 Ratchets A speculum with six ratchets setting for multiple opening options and the extra depth of the ratchets provides extra security to avoid unexpected closing of the speculum whilst on the horse's mouth. Comes with biothene straps which are flexible, easy to clean and use. Available in normal or extra wide bite plates. This product comes with 10 years guarantee.



**Equine Mouth Speculum
Di-0022**

Haumann Speculum This four ratchet speculum is ideal for day to day examination and dentistry procedures. The extra depth of the ratchets locks the speculum securely when horse's mouth is opened. Comes with biothene straps which are water proof and flexible to use. Available in normal or extra wide bite plates. This product comes with 10 years guarantee.



**Equine Mouth Speculum
Di-0023**

Millennium Speculum with Biothene Straps Premium quality speculum has deep but small ratchet settings to ensure a comfortable fit on each individual horse. It can be fitted with normal or extra wide bite plates. Comes with hardwearing biothene straps which are water proof and easy to clean. Available in normal or extra wide bite plates.



**Equine Mouth Speculum
Di-0024**

Pony Speculum with Biothene Straps Ideal for small horses and ponies. This smaller size speculum enables the technicians to work as normal but with more secure fitting



**Equine Mouth Speculum
Di-0025**

Comes with ball and pin option. Insert between cheek teeth to hold incisors apart.



**Equine Mouth Speculum
Di-0026**

Used to extend the equine cheek away From

Equine Mouth Speculums



Equine Mouth Speculum

Di-0027

Equine Dental Head Stand This unit is easy to assemble, dismantle and to store away; it has a unique feature of a spring loaded height adjuster and locking system. The headrest is firmly padded and finished in fine quality leather and additional removable pvc cover (supplied) which is washable to aid hygiene.



Equine Mouth Speculum

Di-0028

Biothene Straps Biothene straps are hard wearing, flexible and easy to use with all types of speculums. Assembled with high quality brass buckles and rivets. Easy to attach to the speculums with the help of screws supplied with the straps.



Equine Mouth Speculum

Di-0029

Incisor Speculum This gag is placed in the mouth of the horse the same way as a bit, It creates easy access to the incisors.

Gritted Floats



Gritted Float
Di-0030

Carbide chip S Float 12" This float comes with a combination of concave and curved ends. They are coated double sided on the curved end and single sided on the concave end. S-Floats come in carbide chip or diamond grit. The body is made from stainless steel and the overall length is 300mm (12") long.



Gritted Float
Di-0031

Diamond Gritted S Float 12" This float comes with a combination of concave and curved ends. They are coated double sided on the curved end and single sided on the concave end. S-Floats come in carbide chip or diamond grit. The body is made from stainless steel and the overall length is 300mm (12") long.



Gritted Float
Di-0032

Handy Float Coated with an fine carbide grit, it has a very sharp cutting profile which creates an excellent finish for general filing and finishing work. The carbide grit form, totally coats both ends of the Handy Float making it extremely easy to control and use in hard to reach places.



Gritted Float
Di-0033

Diamond Gritted S Float 7" Grit diamond plated coated on both side. It has a medium cutting profile to ensure durability and excellent finish for general filing and finishing work.

Scalers & Probes



Scalers & Probes

Di-0034

Dental Explorer Fine Point Overall Length 43cm Fine Dental explorer for examination and exploring of dental teeth.



Scalers & Probes

Di-0035

Dental Scaler Overall Length 43cm Small hook sickle scaler for the removal of unwanted supragingival deposits on teeth. Finer head allows for more precision work or for work in areas hard to reach.



Scalers & Probes

Di-0036

Sickle Scaler Overall Length 43cm A bit larger curved hooksickle scaler for removal of unwanted supragingival deposits from teeth. With super grip handle and manufactured from high quality stainless steel.



Scalers & Probes

Di-0037

Gracey Curvette Left overall Length: 43cm The angled head of the Gracey Curvette allows for better exploration of the horse's mouth. Assists in the removal of unwanted deposits on teeth and for investigative treatments and procedures. With super grip handle and manufactured from high quality stainless steel.



Scalers & Probes

Di-0038


Gracey Curvette Right Overall Length: 43cm The angled head of the Gracey Curvette allows better access to tricky areas within the mouth. Assists in the removal of unwanted deposits on teeth and for investigative treatments and procedures. With super grip handle and manufactured from high quality stainless steel.



Scalers & Probes

Di-0039

Periodontal Probe Banded Overall length: 43cm Banded probe enables the practitioner to judge depths whilst examining the horse's mouth. With super grip handle for expert precision when used. Manufactured from high quality stainless steel.



Scalers & Probes

Di-0040

Periodontal Instrument Set Set includes: 1 peridontal Probe, 1 explorer fine point, 1 sickle scaler large hook, 1 scaler small hook, 1 equine needle pick and 1 dental mirror.